
CE 0197

**Heinen + Löwenstein**  
Arzbacher Strasse 80 56130 Bad Ems Germany  
T: +49-(2603) 960005 F:+49-(2603) 9600-50  
www.hul.de

WM 68110 12/2015 DE, EN, FR, NL, TR, IT, ES(MX), PT(BR), RU, SV, FI






# JOYCEclinic FF

Full Face Mask  
vented, non vented

*Gebrauchsanweisung*  
*Instructions for Use*  
*Mode d'emploi*  
*Gebruiksaanwijzing*  
*Istruzioni d'uso*  
*Kullanma Kılavuzu*

*Manual de instrucciones*  
*Instruções de uso*  
*Инструкция по пользованию*  
*Bruksanvisning*  
*Käyttöohje*

## 1 Bedienung

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## 1 Operation

Use the template on the title page (see Figure 4) of these instructions for use to select the correct mask size. See the following figures for how to put on, adjust and remove the mask:

- 1 Putting on the mask
- 2 Adjusting the mask
- 3 Removing the mask

## 2 Introduction

### 2.1 Intended use

The JOYCEclinic FF full-face mask is used as a connecting element between the patient and the therapy device for non-invasive positive pressure respiratory therapy in hospitals/clinics. It is for a single patient with respiratory insufficiency (> 30 kg) who is breathing spontaneously and who may be subjected to non-invasive ventilation with pressure support. The mask may only be used in combination with therapy devices which have the appropriate alarms and safety systems in case the device fails.

Use JOYCEclinic FF NV only in combination with therapy devices which have an active exhalation valve.

### 2.2 Contraindications

The mask may not be used, or may be used only with particular caution, if any of the following symptoms are present:

erosions and ulcerations, skin allergies, rashes on the face, pressure marks on the face, claustrophobia, anxiety, deformities of the face or nasopharynx, taking of medication which may induce vomiting, the necessity for immediate intubation, pronounced reflux, dysfunctional cough reflex. Please also observe the contraindications in the instructions for use for your therapy device.

### 2.3 Side effects

Nasal congestion, dry nose, dry mouth in the morning, feeling of pressure in the sinuses, irritated conjunctiva, skin rashes, pressure marks on the face, irritating noises when breathing. Masks may cause teeth, gums or jaw to ache.

## 3 Safety

### **Risk of injury from the supply of oxygen!**

Oxygen can become deposited in clothing, bedlinen and hair. In conjunction with smoking, naked flames

or electrical equipment, it can cause fires and explosions.

- ⇒ Do not smoke.
- ⇒ No naked flame.

### **Risk of injury from re-inhaling CO<sub>2</sub>!**

If the mask is used incorrectly, CO<sub>2</sub> may be re-inhaled.

- ⇒ Vented variant only: do not close off the exhalation system of the mask.
- ⇒ Only put on the mask for an extended period if the therapy device is running.
- ⇒ Only use the mask within the quoted therapy pressure range.
- ⇒ Do not use the mask on patients who are unable to take it off themselves.

### **Risk of injury if the mask slips!**

If the mask slips or falls off, the therapy is ineffective.

- ⇒ Monitor patients with restricted spontaneous respiration.
- ⇒ Activate low pressure/leakage alarms on the therapy device.

### **Risk of injury from anesthetic gases!**

Aesthetic gas may escape through the exhalation valve and put third parties at risk.

- ⇒ Never use the mask during anesthesia.

## 4 Product description

### 4.1 Overview

A diagram of the individual parts can be found on the title page.

1. Headgear
2. Forehead cushion
3. Forehead support
4. Retaining ring
5. Elbow
6. Anti-asphyxia valve
7. Mask body
8. Mask cushion
9. Headgear clip
10. Quick-release cord

## 4.2 Compatible devices

With many combinations of device, actual pressure in the mask will not correspond to prescribed therapy pressure, even if the therapy device is displaying the correct pressure. Have the combination of devices adjusted by a physician or specialist dealer so that actual pressure in the mask corresponds to therapy pressure.

## 4.3 Exhalation system

### Vented variants

Vented variants (transparent retaining ring and elbow) have an integrated exhalation system. The retaining ring and mask body are shaped so that there is a gap between these parts. The exhaled air can escape through this gap.

### Non-vented variants

Non-vented variants (blue retaining ring) do not have an exhalation system. Only use non-vented full-face masks in combination with therapy devices which have an active exhalation valve and have alarms and safety systems for any device failure. If you are using the non-vented variant with an external exhalation valve, please see the associated instructions for use.

## 4.4 Anti-asphyxia valve

### ⚠ WARNING

#### Risk of asphyxia if anti-asphyxia valve not working properly!

Residues may cause the valve to stick and lead to CO<sub>2</sub> being re-inhaled.

⇒ Check before every use that the openings of the anti-asphyxia valve are clear.

If the therapy device fails, the anti-asphyxia valve opens so that the patient can breathe ambient air.

## 5 Hygiene treatment

The mask is intended for single use only and should not be disinfected or cleaned after use on the patient. If the mask becomes slightly soiled, we recommend cleaning it with a cloth wetted with isopropyl alcohol (70 % v/v).

### 5.1 Change of patient

The mask is intended for single use only.

## 6 Disposal

Dispose of the mask in accordance with hospital regulations.

## 7 Troubleshooting

Fault	Cause	Remedy
Pain from pressure on the face.	Mask too tight.	Loosen headgear slightly.
Draft in the eye.	Mask too loose.	Tighten headgear slightly.
	Mask does not fit.	Use mask of different size.
Therapy pressure is not reached.	Mask not correctly adjusted.	Re-adjust mask (see Figure 2).
	Patient circuit leaking.	Check connectors, check tubes properly located.

## 8 Technical data

Product class to directive 93/42/EEC	IIa
Dimensions (W x H x D)	
Size S	100 mm x 140 mm x 105 mm
Size M	100 mm x 148 mm x 105 mm
Size L	100 mm x 160 mm x 105 mm
Weight	
Size S	103 g
Size M	105 g
Size L	108 g
Therapy pressure range	0 hPa - 50 hPa
Tube connection: tapered connector to EN ISO 5356-1	Ø 22 mm (female)
Temperature range: operation storage:	+5 °C to +40 °C -20 °C to +70 °C
Flow resistance, vented at 50 l/min at 100 l/min	0.03 hPa 0.1 hPa
Flow resistance, non-vented at 50 l/min at 100 l/min	0.08 hPa 0.32 hPa
Flow resistance, anti-asphyxia valve	
Inspiration at 50 l/min:	0.5 hPa
Exhalation at 50 l/min:	0.5 hPa

Switching pressure	
Anti-asphyxia valve	
• Open:	1 hPa
• Close:	2 hPa
Geometrical dead space:	
Size S	210 ml
Size M	240 ml
Size L	290 ml
Unintended leaks:	
4 - 20 hPa	≤ 4.0 l/min
20 - 50 hPa	≤ 6.0 l/min

No parts of the mask contain latex, PVC (polyvinyl chloride) or DEHP (diethylhexyl phthalate).

# CE 0197

The right to make design modifications is reserved.

## 9 Warranty

The manufacturer Weinmann gives the customer a limited manufacturer warranty on new original Weinmann products, and on any replacement part fitted by Weinmann, in accordance with the warranty conditions applicable to the product in question and in accordance with the warranty periods from date of purchase as listed below. The warranty conditions can be downloaded from [www.weinmann-medical.com](http://www.weinmann-medical.com) on the Internet. We will also send you the warranty conditions on request.

Please bear in mind that any claim to warranty and liability shall be void if neither the accessories recommended in the instructions for use nor genuine replacement parts are used.

In the event of a claim under warranty, please contact Heinen + Löwenstein.

Product	Warranty periods
Masks including accessories	6 months

## 10 Declaration of conformity

The manufacturer Weinmann Geräte für Medizin GmbH + Co. KG (Kronsaalsweg 40, 22525 Hamburg, Germany) hereby declares that the product complies with the relevant regulations of Directive 93/42/EEC governing medical devices. The complete text of the declaration of conformity can be obtained from Heinen + Löwenstein.

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