Allegra M30

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The safety of the **secretion aspirator Allegra M30** complies with the acknowledged rules of technology and meets the requirements of the **German Medical Devices Act**.

The **secretion aspirator Allegra M30** bears the **CE marking CE0044** in accordance with the EU Council Directive 93/42/EEC concerning medical devices and meets the essential requirements of Annex I of this directive.

The **secretion aspirator Allegra M30** has been tested in accordance with IEC 62353.

The **quality management system** applied by Löwenstein Medical GmbH & Co. KG is certified in compliance with the relevant international standards.

The **secretion aspirator Allegra M30** is a medical aspirator classified as class IIa in accordance with EU Council Directive 93/42/EEC, Annex IX.

Errors and omissions excepted.

Contents

I User information.....46

I.1	Using this instructions for use	46
I.2	Icons.....	46
I.2.1	Device, packaging and accessories.....	46
I.2.2	Operating elements.....	47
I.3	Symbol convention.....	47
I.4	Glossary	48
I.5	Purpose.....	49
I.5.1	Intendend use.....	49
I.5.2	Essential features.....	49
I.5.3	Applied parts	49
I.5.4	Indications.....	49
I.5.5	Contraindications.....	49
I.5.6	Restrictions on use	49
I.6	Basic safety instructions - ATTENTION!.....	50
I.7	User requirements.....	52
I.8	Information on product liability.....	52
I.9	Material compatibility	53

2 Product description54

2.1	Whole view of the ME system.....	54
2.2	Scope of delivery, basic device.....	54
2.3	Product properties	54
2.4	Information on the bacterial filter.....	55
2.5	Information on the battery.....	56
2.6	Warranty.....	56

3 Operation.....57

3.1	Set-up and startup.....	58
3.1.1	Connecting the secretion aspirator Allegra M30	58
3.2	Startup.....	59
3.2.1	Connecting the reusable secretion canister.....	59

3.3	Operation.....	62
3.3.1	Control elements and indicators of the secretion aspirator Allegra M30.....	62
3.3.2	Setting the vacuum.....	64
3.3.3	Aspiration.....	64
3.4	Decommissioning.....	65

4 Maintenance.....67

4.1	Cleaning and care in outpatient and inpatient care.....	67
4.1.1	General information.....	67
4.1.2	Cleaning and disinfection of the surface of the device.....	67
4.1.3	Cleaning and disinfection of the reusable secretion canister.....	68
4.1.4	Cleaning/disposal of the suction tube.....	68
4.1.5	Cleaning and disinfection of the rinsing bottle.....	68
4.2	Cleaning and care in the homecare area.....	69
4.2.1	General information.....	69
4.2.2	Cleaning of the surface of the device.....	69
4.2.3	Hygienic cleaning and disinfection of the reusable secretion canister.....	70
4.2.4	Hygienic cleaning and disinfection of the suction tube.....	70
4.2.5	Hygienic cleaning and disinfection of the rinsing bottle.....	71
4.3	Further use of the device.....	71
4.4	Maintenance and service.....	71
4.5	Inspection of the secretion aspirator Allegra M30.....	72

5 Troubleshooting73

6 Transport, storage and disposal.....74

6.1	Decontamination prior to shipment.....	74
6.2	Storage.....	74
6.3	Disposal.....	74

7 Technical data75

8 EMC information.....76

8.1	Electromagnetic environment in which the secretion aspirator Allegra M30 may be operated.....	76
8.2	Handling of electromagnetic interactions.....	79
8.3	List of all cables and transformers replaceable by the user.....	79

9 Order information accessories80

10 Publishing information.....81

I User information

I.1 Using this instructions for use

Please read this entire instructions for use before operating the **secretion aspirator Allegra M30** for the first time. Alongside instructed medical professionals and trained relatives of patients, it is also intended that patients operate the device themselves, providing that they have also been instructed accordingly. In this case all functions can be used safely. If support is required for startup, operation or maintenance, please contact Löwenstein Medical GmbH & Co. KG (see section I 0). Please also report unexpected operating performance or incidents to Löwenstein Medical GmbH & Co. KG.

Please read the safety instructions (see section I.6) to avoid hazards.

This instructions for use forms an integral part of the **secretion aspirator Allegra M30**. Therefore please keep this instructions for use in an easily accessible place.

Include this instructions for use when passing the **secretion aspirator Allegra M30** on to third parties.

I.2 Icons

I.2.1 Device, packaging and accessories

Symbol	Meaning
	Attention: possible bodily injury, health risks or possible property damage
	Follow the instructions for use
	Power supply unit
	Protection class II
	Protection class: Type BF (Body Floating)
	Temperature limitation
	Air pressure limitation
	Humidity limitation
	Do not use if packaging is damaged!
	International Protection: IP22 (see section I.4)
	Note with useful information and tips

Symbol	Meaning
	This device must not be disposed of in domestic waste.
	Order number
	Serial number
	Lot number
	Date of manufacture
	Manufacturer
	Operating time (power-on time and power-off time)
	Do not reuse.

1.2.2 Operating elements

Symbol	Meaning
	Vacuum setting "Low" (≤ -0.2 bar)
	Vacuum setting "Medium" (≤ -0.5 bar)
	Vacuum setting "High" (≤ -0.9 bar)
	LED status indicator for battery operation
	LED status indicator for mains operation
	On/off button

1.3 Symbol convention

Symbol	Meaning
•	Enumeration
1. 2.	Perform the actions in the specified order.

I.4 Glossary

A	Aspirate	Aspirate is the generic term for secretions, bodily fluids and flushing liquids that are typically accumulated during aspiration of the upper airways. These can be aspirated easily using the device described here.
C	Contamination	Contamination means that bacteria and viruses from the aspirate have come into contact with the device.
D	Degree of protection	The degree of protection indicates the protection of applied parts against electric shock. Applied parts of type BF have to be installed isolated from earth and are not suitable for direct application to the heart.
I	incl.	Abbreviation for "including"
	IP22	International Protection/Protection Class The Protection Class specifies the degree of protection of the device against contact and ingress of liquids. The secretion aspirator Allegra M30 is protected against finger access and falling water drops at an inclination of up to 15°.

M	ME	Medical electrical
	MRI	Abbreviation for "magnetic resonance imaging" This technique generates sectional images of the human body with the aid of a very strong magnetic field for the purpose of analyzing the organs.
O	Overflow	Overflow means that aspirate is sucked into the interior of the device.
P	Processing	The processing procedure is required for each new patient. The term processing means that all parts that have or may have come into contact with the aspirate must be cleaned, disinfected and replaced if necessary. Reconditioning may only be performed by Löwenstein Medical GmbH & Co. KG or by a service partner authorized by Löwenstein Medical GmbH & Co. KG.
R	resp.	Abbreviation for "respectively"

1.5 Purpose

1.5.1 Intended use

The **secretion aspirator Allegra M30** is a network-independent mobile medical aspirator and serves for temporary and primarily spontaneous aspiration of aspirate from the trachea.

Typical areas of application are:

- in outpatient and inpatient care (professional healthcare facility environment) and in the homecare area (home healthcare environment).

1.5.2 Essential features

- Generation of vacuum (high vacuum)
- Generation of flow (high flow)

1.5.3 Applied parts

There are no applied parts, however housing and suction tube are considered as applied parts of type BF.

1.5.4 Indications

- Aspiration on impaired respiratory function
- Aspiration on blood, secretion and food residues from the oral cavity, the pharyngeal zone and the bronchial system

1.5.5 Contraindications

The **secretion aspirator Allegra M30** is contraindicated for the following applications:

- Liposuction
- Gynecology applications
- Dental applications
- Thoracic drainage
- Continuous drainage
- Applications in wound areas

1.5.6 Restrictions on use

- In medical rooms where potential equalization is required (e.g. heart surgery)
- In hazardous areas/in the MRI environment
- Outdoors/during transport

1.6 Basic safety instructions - ATTENTION!

Warning of damage due to improper power supply

Improper operation causes overvoltage in the device, which may be transmitted to the operator:

- Ensure prior to startup that the mains supply for connecting the **secretion aspirator Allegra M30** is suitable at 100 V to 240 V AC and a mains frequency of 50-60 Hz.
- Ensure prior to startup in UL listed markets such as the USA and Canada that the mains supply is designed for a supply voltage of 120 V AC.
- Only use the supplied power supply unit with the **secretion aspirator Allegra M30** (type: GTM91099-6015-3.0-T2 from GlobTek Inc.).

Risk of application under conditions that are not approved

- The device is not intended to be used in medical rooms where potential equalization is required (e.g. heart surgery).
- The device is not intended to be used in hazardous areas/in the MRI environment.
- The device is intended for mobile, temporary outdoor use (in accordance with IP protection class 22), but not for permanent outdoor location and operation.

Health risks due to the handling of infectious or pathogenic germs

Infectious and pathogenic germs in the aspirate cause health risks.

- Never operate the device without bacterial filter.
- Prior to each application, make sure that the bacterial filter is dry and clean to ensure faultless operation. Pay attention that the bacterial filter does not get moist during the performance of cleaning/disinfection.
- Always perform the aspiration using a suitable sterile disposable catheter. The suction tube must never come into direct contact with the aspiration site.
- Follow the hygiene, cleaning and decontamination instructions.

Hazard of persons due to strangulation

- People may strangle themselves on the tubing or the mains cable, especially if tubes resp. cables are unduly long.
- Ensure that no unauthorized/uninvolved persons are near the device during aspiration.
- Store the device incl. accessories in the shipping carton until the next time of use.

Risk of damage due to electromagnetic phenomena

Medical electrical devices are subject to special precautionary measures regarding electromagnetic compatibility and must be installed and operated in accordance with the EMC information provided in the accompanying documents (see section 8).

Hazard of persons due to improper handling

- Only use the device for its intended purpose.
- Never use the device for aspiration in the low vacuum range (e.g. thoracic drainage).
- If aspiration procedures are carried out too often, slight bleeding may occur.
- When using the device for the treatment of children, the vacuum setting "Low" must be used.
- When using the power supply unit, ensure that it is connected to the aspirator firstly and afterwards to the mains supply (100 V to 240 V AC). The disconnection of the power supply unit from the mains supply must be performed in exactly the opposite sequence (first remove the power supply unit from the mains supply (100 V to 240 V AC) and then remove it from the aspirator).
- Never touch parts of other devices than ME devices in the patient environment and the patient simultaneously.

Safety defects due to improper accessories and spare parts

The use of accessories and spare parts, detachable parts or materials that are not recommended by Löwenstein Medical GmbH & Co. KG and specified in the instructions for use may compromise the safety and function of the device. Any warranty is excluded for damage caused by using accessories and spare parts not recommended or by improper use.

- Only use original accessories and spare parts.

Warning of safety defects due to improper connections of the ME system

The connection of the ME system with other devices or installations resp. pieces of equipment not recommended by Löwenstein Medical GmbH & Co. KG and not specified in the instructions for use may compromise the safety and function of the ME system. Any warranty is excluded for damage caused by connecting devices or installations resp. pieces of equipment not recommended with the **secretion aspirator Allegra M30** or by improper use.

- Only connect recommended original parts with the **secretion aspirator Allegra M30**.

Damage of the device by heat

- Do not cover the power supply unit.
- Keep the aspirator as well as the mains cable and the power supply unit away from other heat sources.
- Always position the aspirator a sufficient distance away from other devices which generate heat (e.g. not directly in front of a cooling fan opening) to prevent the aspirator heating up excessively.

Damage of the device due to improper handling

- Never aspirate flammable, corrosive or explosive liquids or gases.
- Do not drop the device.
- Prior to each use check the housing for any damage and do not operate the device in case of apparent housing damage.
- Prior to each use check the suction tube, secretion canister and bacterial filter and if applicable any other accessories, that are subject to wear and damage, to ensure that the components are in a perfect condition and proper operation of the device is guaranteed. If this is not the case, replace the parts immediately.

Inspection of the internal power supply

Since the internal battery is not maintained automatically in a fully operational state, it is necessary to check the charging status regularly and if necessary, the battery must be replaced by service staff. The battery may only be replaced by trained service staff, since the replacement by inadequately trained personnel could result in a hazard (such as excessive temperatures, fire or explosion)!

Possible physiological effects and unobvious risks

- Select the vacuum range depending on the respective patient and on the medical indication to avoid personal injuries.
- During the aspiration procedure, ensure that the device is in an upright position to prevent the secretion canister overflowing.
- Always place the device upright on a sturdy and flat base without slope. Make sure that the device cannot be knocked over or fall in such a way that persons could be hit by the falling device.
- Only operate the device with inserted vacuum pump and closed lid. The lid above the vacuum pump may only be opened by adequately trained service staff! Likewise, the vacuum pump may only be replaced by adequately trained service staff.
- Other devices, examinations or treatments may possibly be influenced by the device. For this reason, special attention should always be paid to other devices as well as to examinations and treatments performed in parallel so that any influence is detected as quickly as possible.
- When using the device, adequate lighting must always be ensured so that all labels can be identified clearly.
- Small, detached parts could be inhaled or ingested. Therefore ensure that no unauthorized persons, children or pets are near the device.
- Although the materials used have been tested for compatibility, in exceptional cases allergic reactions to the exposed materials on the device and its accessories may occur. This applies especially for contact injuries after prolonged contact. In such cases, turn to a doctor immediately.

Damage of the device by ingress of liquids

The **secretion aspirator Allegra M30** has the IP classification IP22 against the ingress of liquids. Nevertheless, protect the device from moisture.

- Do not use the device near splashing water.
- Do not use the device in damp rooms where humidity is >90 % (no condensation), in the bath or when showering (risk of electric shock).
- Avoid wet conditions on power supply unit, control panel and socket for the power supply unit.
- Never submerge the device in water or other liquids (also not while it is switched off).

Known, identifiable or foreseeable conditions for medical care within a domestic environment

- Children and pets should be kept away from the device. Make sure that the device cannot be knocked over or fall in such a way that persons could be hit by the falling device.
- Prior to connecting the power supply unit, ensure that the voltage of the device corresponds to the domestic power supply.
- The device may not be used in damp rooms where humidity is >90 % (no condensation). Avoid wet conditions on power supply unit, control panel and socket for the power supply unit.
- Ensure that fluff and dust are immediately removed from the device incl. accessories, in order to ensure unimpaired functionality. Furthermore no pests may be allowed in the proximity of the device, as there is a risk they may get inside the device causing damage.
- Never expose the device incl. accessories to direct sunlight as this could result in a strong heating effect, causing functional impairments.
- Some devices and sources that are commonly used in the domestic environment may cause potential faults in the device incl. accessories, e.g. fireplaces or radiant heaters (strong heating of the device) resp. inhalers or steam kettles (excessive air humidity). Do not operate such devices and sources in the vicinity of the **secretion aspirator Allegra M30**.

I.7 User requirements

The **secretion aspirator Allegra M30** may only be operated and used by instructed and appropriately trained persons. Familiarize yourself with the mode of operation of the **secretion aspirator Allegra M30** prior to startup.

Training on how to handle the **secretion aspirator Allegra M30** is provided by Löwenstein Medical GmbH & Co. KG or an authorized distribution partner of Löwenstein Medical GmbH & Co. KG. Product training, among other things, includes an explanation of the design and function of the device, the handling of the device, the cleaning and disinfection as well as the procedure to be followed for each new patient and for disposal.

I.8 Information on product liability

The liability for the operation of the device is channeled to the operator in the following cases:

- the device is used outside its intended use,
- the device is not used in accordance with the instructions for use,
- the device is opened by unauthorized personnel,
- the safety seal is removed/broken,
- installation, settings, enhancements, routine maintenance or repairs are performed by unauthorized personnel,
- original accessories and spare parts have not been used.

Advice for the responsible organization:

The assembly of ME systems and modifications during their actual service life shall be evaluated based on the requirements of the applying standards.

1.9 Material compatibility


Aggressive substances can damage the device and the accessories.

Follow the instructions for cleaning and care (see section 4.1 resp. 4.2).

2 Product description

2.1 Whole view of the ME system

- A Secretion canister (reusable) with integrated overflow protection in the lid (fig. 19)
B Control panel
C On/off button
D Socket for power supply unit (12V DC)
E Power supply unit (GlobTek GTM91099-6015-3.0-T2)
F Mains cable
G Suction tube
H Lid vacuum pump
I Carrying handle

2.2 Scope of delivery, basic device

- the **secretion aspirator Allegra M30**
- the instructions for use of the **secretion aspirator Allegra M30**
- the power supply unit GlobTek GTM91099-6015-3.0-T2
- reusable secretion canister system (incl. suction tube and fingertip)
- rinsing bottle, reusable (250 ml)
- charging instructions, multilingual
- safety instruction for handling battery packs
- "Used Medical Device" label and decontamination certificate
- test report in accordance with IEC 62353

2.3 Product properties

Hazards of persons due to improper handling

- Only use the device for its intended purpose.
- Never use the device for aspiration in the low vacuum range (e.g. thoracic drainage).

Damage to the device due to improper handling

- Never aspirate flammable, corrosive or explosive liquids or gases.
- Do not drop the device.
- Do not operate the device in case of apparent housing damage.

The **secretion aspirator Allegra M30** is an especially handy, battery-operated medical aspiration device for mobile and stationary use.

The **secretion aspirator Allegra M30** has a flow rate of max. 30 l/min $\pm$ 2 l/min (see section 7 "Technical data").

The **secretion aspirator Allegra M30** is operated via the internal battery or via the supplied power supply unit that can also be used to recharge the battery. A rotary vane pump generates the desired vacuum.

After 5 minutes of operation, the device switches off automatically to prevent overheating. The device should be left switched off for at least 5 minutes to ensure sufficient cooling of the device. In emergency situations, however, it is possible to switch the device back on immediately several times.

Moreover an overheating protection prevents the battery from overheating by interrupting the charging process if the temperature of the battery is too high (e.g. due to unfavorable environmental conditions).

After it is switched on, the vacuum pump generates a vacuum in the secretion canister, which is used to aspirate secretion (via a suitable sterile aspiration set). The aspirate is directed away from the patient and collected in the secretion canister. If the secretion canister is full, the aspiration is stopped by means of a mechanical overflow protection integrated into the lid of the secretion canister.

The desired vacuum is set by buttons in the control panel (fig. 1 (B), see section 1.2.2)). The secretion canister is attached by sliding it directly into the housing structure of the aspirator (fig. 8, 9). This also establishes the vacuum connection between the aspirator and secretion canister.

The device is charged or operated via the supplied 12V power supply unit (fig. 1 (E)).

2.4 Information on the bacterial filter

Health risks due to the handling of infectious or pathogenic germs

Infectious and pathogenic microbes in the suction material can pose a risk to health.

- Never operate the device without bacterial filter.
- Prior to each application, make sure that the bacterial filter is dry and clean to ensure faultless operation.

The bacterial filter integrated between the vacuum pump and the secretion canister consists of a hydrophilic (self-sealing) material and has an efficiency of 99 % against bacteria and viruses. Thus it effectively protects the vacuum pump against overflow and contamination and enables fast, simple and cost-effective processing.

The **integrated bacterial filter** is only intended for **use on one patient**. It must therefore be replaced by the operator each time a new patient uses the device (see section 3.2.1). If the device is used on one patient exclusively, the bacterial filter must be replaced at least every 2 weeks.

To ensure consistent performance, the integrated bacterial filter must also be replaced by authorized service staff if it has come into **contact with the aspirate** (blocked) or **during maintenance/repair**.

Since the drive unit is mechanically isolated from the vacuum pump, it cannot become contaminated during use independent from the bacterial filter.

2.5 Information on the battery

The charge state of the battery is indicated by the corresponding LED status indicator (for explanation see section 3.3.1).

Prior to the initial startup of the **secretion aspirator Allegra M30** it is urgently recommended to charge the battery completely and to repeat this after the first uses.

The **secretion aspirator Allegra M30** is equipped with a lithium ion battery which only has a limited self-discharge on its disposal in contrast to conventional battery types.

Lithium ion batteries do not have a memory effect. Therefore they can and should be recharged at any time after initial charging. Only frequent short charging should be avoided. The typical service life of the battery amounts to approximately 500 load cycles, subsequently approximately 70% of the primary capacity are still available.

The battery of the **secretion aspirator Allegra M30** is protected against depth discharge by means of protective measures; however the instructions on charging mentioned above should be observed. The battery is also protected against overheating during charging. If the battery temperature is exceeded during charging due to ambient conditions not intended, charging will be automatically interrupted temporary to allow cooling. This feature ensures the safety of the battery and protects it from wear.

2.6 Warranty

WARNING:

Modifications of the ME device are not permitted.

The **secretion aspirator Allegra M30** is covered by warranty for 2 years. No extension or renewal of the warranty is granted in the event that warranty work is performed.

The battery is covered by warranty for 6 months. Wearing parts are excluded from the warranty.

Löwenstein Medical GmbH & Co. KG is responsible for impacts on safety, reliability and specific performance only if:

- original accessories and replacement parts are used,
- maintenance and repair are performed by professionals authorized by Löwenstein Medical GmbH & Co. KG or by Löwenstein Medical GmbH & Co. KG itself,
- the **secretion aspirator Allegra M30** is used and operated in accordance with the instructions for use and within its intended use.

All warranty claims shall be void if:

- the device is opened by unauthorized personnel,
- the safety seal is removed/broken,
- repairs are performed by unauthorized personnel,
- modifications are made to the device,

since the basic safety of the device can no longer be guaranteed in these cases and functional limitations may occur.

3 Operation

Hazard of persons due to improper handling

- Please read through sections 3.1 and 3.2 carefully!
- Only perform aspiration in the respiratory tract after instruction by specialized personnel!
- Only use suitable sterile disposable catheters for aspiration! Commercially available sterile disposable catheters from CH10 to CH16 that are intended for the use in the tracheal area can be connected.
- If aspiration procedures are carried out too often, slight bleeding may occur.
- When using the device for the treatment of children, the vacuum setting "Low" must be used.

Malfunction due to aspirated secretion

- Change the reusable secretion canister when it is half-full to curtail frothing. If the reusable secretion canister is full, the mechanical overflow protection is triggered. This disrupts the aspiration process.
- Switch off the device when emptying the reusable secretion canister.
- If secretion has been aspirated into the device, it must be properly processed by Löwenstein Medical GmbH & Co. KG or by a service partner authorized by Löwenstein Medical GmbH & Co. KG!

Damage of the device due to insufficient acclimatization

- After the device has been exposed to temperatures between -25°C and +70°C during transport/storage according to the technical data (see section 7), it must firstly acclimatize for ca. 2 h at room temperature (ca. 20°C) before the intended use is possible.

3.1 Set-up and startup

The following section describes the operating elements, connections and the startup of the **secretion aspirator Allegra M30**:

3.1.1 Connecting the secretion aspirator Allegra M30

Known or identifiable conditions for medical care within a domestic environment

- Children and pets should be kept away from the device to ensure that it cannot be knocked over or dropped.
- Prior to connecting the power supply unit, ensure that the voltage of the device corresponds to the domestic power supply.
- Do not use the device in damp rooms where humidity is >90 % (no condensation), in the bath or when showering (risk of electric shock). Avoid wet conditions on power supply unit, control panel and socket for the power supply unit. Never submerge the device in water or other liquids (also not while it is switched off).

Check the power supply unit incl. mains cable for potential damages before each use and replace it if there is any damage.

Use the socket for the 12V power supply unit (fig. 1(D)) to connect the device to the mains power supply with the aid of the supplied power supply unit to charge or operate it as needed. Make sure that the device is positioned in such a way that it is possible to easily disconnect it later.

Information on the permissible environmental conditions during operation can be found in section 7 "Technical data".


Fig. 2


Connection of the secretion aspirator Allegra M30 to patient and accessories

3.2 Startup

- Take the device and the accessories out of the packaging.
- During the aspiration procedure, ensure that the device is in an upright position to prevent the secretion canister overflowing.
- Prior to the initial startup it is essential to follow the safety instructions in section 1.6.
- Only operate the device with inserted bacterial filter and connected secretion canister. Bacterial filter and secretion canister are already fitted on delivery. Note the following instructions for all further applications.
- Always have one backup bacterial filter at hand, since it is absolutely essential for safe operation!
- Prior to each application, check that the reusable secretion canister is completely clean to avoid frothing.

3.2.1 Connecting the reusable secretion canister

1. Fold the moving part of the secretion canister lid slightly upwards and place the secretion canister lid on the secretion canister.



3 Operation

2. Fold the moving part of the secretion canister lid downwards so that the hooks snap into place and the secretion canister is securely closed.



3. Insert the bacterial filter into the device by introducing the filter with the lever pointing vertically downwards in the opening provided. Subsequently move the lever a quarter turn clockwise so that the lever is in the horizontal position afterwards.


Please note that the reusable secretion canister cannot be connected if the bacterial filter has not been inserted correctly.

4. Insert the secretion canister into the canister fitting provided on the device by guiding it over the rack on the bottom and pivoting it into the device until it snaps into place.


5. Connect the suction tube to the port in the middle of the secretion canister lid. In doing so, ensure that the tube is fixed sufficiently tight.

6. Connect the other ending of the suction tube to the aspiration catheter with the aid of a fingertip.

3.3 Operation

3.3.1 Control elements and indicators of the secretion aspirator Allegra M30

Symbol	Description	Meaning of display
	button illuminates continuously green	device is running, vacuum setting "High" is adjusted (≤ -0.9 bar)
	button illuminates continuously red	canister full/filter blocked (motor switches off automatically)
	button illuminates continuously green	device is running, vacuum setting "Medium" is adjusted (≤ -0.5 bar)
	button illuminates continuously red	canister full/filter blocked (motor switches off automatically)
	button illuminates continuously green	device is running, vacuum setting "Low" is adjusted (≤ -0.2 bar)
	button illuminates continuously red	canister full/filter blocked (motor switches off automatically)
	display illuminates continuously green	battery capacity 100% to 60%
	display illuminates continuously orange	battery capacity $< 60\%$ to 10%
	display illuminates continuously red	Battery capacity $< 10\%$, connect power supply unit immediately. Switch off the device after completion of the aspiration procedure for charging!
	display blinks green	Battery is being charged.
	display blinks red	Overttemperature protection, interruption of charging process until battery temperature returns to permissible range.
		Battery defective, operation only possible with connected power supply unit. Please contact your service partner!
	display illuminates continuously green	Power supply unit is connected, voltage available.
	display blinks red	Incorrect power supply unit is connected, no charging of battery possible!
	button illuminates continuously green	Device is switched on and ready for use.
	button illuminates continuously orange	Device still runs for 30 sec. before it switches off automatically. Afterwards the button will still be illuminated for 10 seconds.
	button illuminates continuously red	Internal error; operation is not possible. Please contact your service partner!


3.3.2 Setting the vacuum

1. Switch on the device by pressing the on/off button (fig. 11(F)).
2. Select the desired vacuum range by pressing the corresponding button ("Low", "Medium" or "High") on the control panel (fig. 11(A, B, C)).

As described in section 3.3.1, the vacuum setting buttons correspond to certain vacuum ranges.

3. You can perform fine adjustments of the set vacuum via the additional air opening on the fingertip (fig. 12).

4. Switch off the device by holding the on/off button (fig. 11(F)) pressed for 1-2 seconds.

3.3.3 Aspiration

1. Perform the aspiration procedure as you were instructed by specialized personnel.
2. Control the aspiration procedure and the desired vacuum via the additional air opening on the fingertip (fig. 12).
3. Decommission the device after the aspiration procedure according to section 3.4.

In any case, please take note of section 2.4 and the important information on the handling of the bacterial filter! Ensure that the battery has a sufficient capacity during the aspiration procedure. The battery status display must signalize green or orange. If the display turns red, the power supply unit must be plugged in immediately. Only in this case the aspiration can be continued without any interruption.

3.4 Decommissioning

Health risks due to the handling of infectious or pathogenic germs



Infectious and pathogenic germs in the aspirate cause health risks.

- Replace the bacterial filter at least every 2 weeks when used exclusively for one patient!
- Wear suitable disposable gloves in doing so.
- The reuse of bacterial filters for multiple patients is prohibited for hygiene and safety reasons!
- For each new patient, processing by Löwenstein Medical GmbH & Co. KG or by an authorized service partner of Löwenstein Medical GmbH & Co. KG is strictly required!
- After each aspiration, all components that have come into contact with the aspirate must be cleaned, disinfected or disposed of.
- The aspirate and the parts contaminated with aspirate must be disposed of properly.

1. Switch off the device after the aspiration by holding the on/off button (fig. 11 (F)) pressed for 1-2 seconds.
2. Disconnect the power supply unit from the mains supply (100 V to 240 V AC) if attached and remove the device plug from the aspiration device afterwards.
3. Remove the suction tube (fig. 1 (G)) from the aspiration catheter and from the secretion canister.
4. Clean or dispose of the suction tube properly (see section 4.1.4 resp. 4.2.4).
5. Remove the reusable secretion canister from the device by pressing the unlocking mechanism on the secretion canister lid and pivoting the secretion canister out of the device afterwards.



6. Fold the moving part of the secretion canister lid upwards so that the hooks are released and take off the lid from the secretion canister.


7. Empty and clean the reusable secretion canister and the lid according to section 4.1.3 resp. 4.2.3.

8. Remove the used bacterial filter before a new patient uses the device or when it has already been used for one patient for 2 weeks. Move the lever that is in the horizontal position a quarter turn counterclockwise so that the lever points vertically downwards subsequently.
9. When the lever is in this position, pull the bacterial filter out of the fitting and dispose of the bacterial filter properly.

10. Insert a new bacterial filter (see section 3.2.1).

11. Clean the surface of the device according to section 4.1.2 resp. 4.2.2.

Recommendation in the event of an extended storage period: Keep the device in suitable packaging, e.g. the shipping carton, until it is used the next time.

4 Maintenance

4.1 Cleaning and care in outpatient and inpatient care

4.1.1 General information

Health risks due to the handling of infectious or pathogenic germs

Infectious and pathogenic germs in the aspirate cause health risks.

- Never operate the device without bacterial filter.
- Replace the bacterial filter at least every 2 weeks when used exclusively for one patient!
- Wear suitable disposable gloves in doing so.
- The reuse of bacterial filters for multiple patients is prohibited for hygiene and safety reasons!
- For each new patient, processing by Löwenstein Medical GmbH & Co. KG or by an authorized service partner of Löwenstein Medical GmbH & Co. KG is strictly required!
- After each aspiration, all components that have come into contact with the aspirate must be cleaned, disinfected or disposed of.
- The aspirate and the parts contaminated with aspirate must be disposed of properly.

Health risks due to the handling of disinfectants

- The use of appropriate protective clothing during disinfection is recommended.
- Follow the disinfectant manufacturer's instructions.

Damage of the device due to improper cleaning agents

- Do not use disinfectants containing acetone. These may damage or disfigure the housing parts as well as the secretion canister lid.
- Follow the instructions for use provided by the manufacturers of the used disinfectants, particularly concerning the compatibility with materials and surfaces as well as the concentration information.
- Löwenstein Medical GmbH & Co. KG recommends "Sekusept® aktiv" for immersion disinfection of the accessories and "Incidin® Plus" or "Incidin® Liquid" for wipe disinfection of the device.

4.1.2 Cleaning and disinfection of the surface of the device

Clean the surfaces of the device regularly and disinfect them at least weekly.

- The device can be wiped with a damp, lint-free cloth.
- Follow the instructions in section 4.1.1 for wipe disinfection.

Repeated cleaning and disinfection procedures may result in minor discolorations of the plastic components of the housing. However these do not affect the function of the device.

If the interior of the device comes into direct contact with liquids or solids, the device must be inspected by Löwenstein Medical GmbH & Co. KG or by a service partner authorized by Löwenstein Medical GmbH & Co. KG.

4.1.3 Cleaning and disinfection of the reusable secretion canister


Fig. 19

- A Reusable secretion canister lid
B Overflow protection with sphere

1. Empty the reusable secretion canister and dispose of the aspirate properly.
2. Remove the sphere before performing the following cleaning instructions (fig. 19(B)).
3. Rinse all components of the reusable secretion canister under running water.
4. Immerse all components of the reusable secretion canister in disinfectant solution, considering the specific concentration quantity.
5. Afterwards, rinse the components thoroughly and let them dry.
6. The reusable secretion canister can also be autoclaved for 20 minutes at 121°C, in that case please remove the overflow protection in advance (fig. 19(B)).
7. After autoclaving, fix the overflow protection in the lid considering the orientation of the opening (fig. 19(B)).
8. Reinsert the sphere according to fig. 19.

Löwenstein Medical GmbH & Co. KG recommends to replace the canister system **after 12 months** if it is used and disinfected frequently!

4.1.4 Cleaning/disposal of the suction tube

1. Dispose of the suction tube properly before the device is used by a new patient!
2. If the device is used for one patient only, rinse the suction tube with clear water (max. 0.5 l) using the rinsing bottle after each aspiration and place it daily in the disinfectant solution recommended by Löwenstein Medical GmbH & Co. KG for immersion disinfection (see section 4.1.1).
3. Afterwards, rinse the suction tube thoroughly with clear water and let it dry.

Even when the suction tube is used for one patient only, Löwenstein Medical GmbH & Co. KG recommends to replace the suction tube at least **every 4 weeks** since the material may become discolored and brittle due to frequent applications!

4.1.5 Cleaning and disinfection of the rinsing bottle

1. Remove the lid of the rinsing bottle and rinse both parts under running water.
2. Immerse the rinsing bottle and the lid in disinfectant solution, considering the specific concentration quantity.
3. Afterwards, rinse the components thoroughly and let them dry.
4. Place the lid back on the rinsing bottle.

Löwenstein Medical GmbH & Co. KG recommends to replace the rinsing bottle **every 4 weeks** if it is used and disinfected frequently!

4.2 Cleaning and care in the homecare area

4.2.1 General information

Health risks due to the handling of infectious or pathogenic germs

Infectious and pathogenic germs in the aspirate cause health risks.

- Never operate the device without bacterial filter.
- Replace the bacterial filter at least every 2 weeks when used exclusively for one patient!
- Wear suitable disposable gloves in doing so.
- The reuse of bacterial filters for multiple patients is prohibited for hygiene and safety reasons!
- Prior to passing the device on to another patient resp. to other persons, processing by Löwenstein Medical GmbH & Co. KG or by an authorized service partner of Löwenstein Medical GmbH & Co. KG is strictly required! For that purpose observe the information regarding the further use of the device (see section 4.3).
- After each aspiration, all components that have come into contact with the aspirate must be cleaned, disinfected or disposed of.
- The aspirate and the parts contaminated with aspirate must be disposed of properly.

Health risks resp. damages of the device due to the handling of disinfectants

- If you have the possibility to use disinfectants, please observe the information in section 4.1. Otherwise please follow the indications of section 4.2.2 till 4.2.5!

4.2.2 Cleaning of the surface of the device

Clean the surfaces of the device regularly.

The device can be wiped with a damp, lint-free cloth. Repeated cleaning procedures may result in minor discolorations of the plastic components of the housing. However these do not affect the function of the device.

If the interior of the device comes into direct contact with liquids or solids, the device must be inspected by Löwenstein Medical GmbH & Co. KG or by a service partner authorized by Löwenstein Medical GmbH & Co. KG.

4.2.3 Hygienic cleaning and disinfection of the reusable secretion canister


Fig. 20

- A Reusable secretion canister lid
- B Overflow protection with sphere

1. Empty the reusable secretion canister and dispose of the aspirate properly.
2. Remove the sphere before performing the following cleaning instructions (fig. 20(B)).
3. Rinse all components of the reusable secretion canister under running water.
4. Place all components of the reusable secretion canister in a water bath at 65°C temperature for 15 minutes once a day.
5. Afterwards, rinse the components thoroughly and let them dry.
6. Reinsert the sphere according to fig. 20.

Löwenstein Medical GmbH & Co. KG recommends to replace the canister system **after 12 months** if it is used and cleaned/disinfected hygienically frequently!

4.2.4 Hygienic cleaning and disinfection of the suction tube

1. Dispose of the suction tube properly before the device is used by a new patient!
2. If the device is used for one patient only, rinse the suction tube with clear water (max. 0.5 l) using the rinsing bottle after each aspiration and place it in a water bath at 65°C temperature for 15 minutes once a day.
3. Afterwards, rinse the suction tube thoroughly with clear water and let it dry.

Even when the suction tube is used for one patient only, Löwenstein Medical GmbH & Co. KG recommends to replace the suction tube at least **every 4 weeks** since the material may become discolored and brittle due to frequent applications!

4.2.5 Hygienic cleaning and disinfection of the rinsing bottle

1. Remove the lid of the rinsing bottle and rinse both parts under running water.
2. Place the rinsing bottle and the lid in a water bath at 65°C temperature for 15 minutes once a day.
3. Afterwards, rinse the components thoroughly and let them dry.
4. Place the lid back on the rinsing bottle.

Löwenstein Medical GmbH & Co. KG recommends to replace the rinsing bottle **every 4 weeks** if it is used and disinfected frequently!

4.3 Further use of the device

The **secretion aspirator Allegra M30** is suitable for further use. However, prior to passing the device on to other patients resp. persons, it has to be processed professionally. For this purpose, please hand the **secretion aspirator Allegra M30** over to Löwenstein Medical GmbH & Co. KG or to qualified personnel authorized by Löwenstein Medical GmbH & Co. KG. In this regard, please observe the indications in section 6.1!

4.4 Maintenance and service

Risk of impairments of performance due to aging effects

Only use the device within its indicated service life of 5 years. If you notice impaired performance or other problems, follow the instructions on troubleshooting in section 5 and contact the service if necessary.

Risk due to maintenance and service during operation

Maintenance and service must not be performed while the device is in use! The device must be switched off before maintenance and service work are performed.

Risk due to the performance of maintenance and service by unauthorized persons

Maintenance and service may only be performed by specialized personnel authorized by Löwenstein Medical GmbH & Co. KG or by Löwenstein Medical GmbH & Co. KG itself.

Risk due to the performance of maintenance and service without sufficient documents

Maintenance and service may only be performed considering the service instructions for the **secretion aspirator Allegra M30**. Löwenstein Medical GmbH & Co. KG will also provide additional documents on request, in case these can support the service staff.

WARNING: Modifications of the ME device are not permitted.

- The **secretion aspirator Allegra M30** is maintenance-free if used in accordance with the instructions for use, exempt from components with limited durability.
- Perform a visual and functional check prior to each application. Also include the accessories of the device in doing so.
- The device may only be opened and repaired by Löwenstein Medical GmbH & Co. KG or by specialized personnel authorized by Löwenstein Medical GmbH & Co. KG in compliance with the service documents specified by the manufacturer as well as with technical and hygienic safety measures.
- The device may be returned to Löwenstein Medical GmbH & Co. KG for repair either directly or via the specialty dealer you purchased the device from.
- All accessories must be cleaned and disinfected before the device is returned. The device itself may be treated with a surface disinfectant.
- In order to avoid time lags, notify Löwenstein Medical GmbH & Co. KG prior to returning the device and mark the shipping carton of the device with the label "Used Medical Device".

4.5 Inspection of the secretion aspirator Allegra M30

Löwenstein Medical GmbH & Co. KG offers its partners and customers fast and professional processing as well as required testing services.

Löwenstein Medical GmbH & Co. KG neither ensures proper operation for the **secretion aspirator Allegra M30** nor is Löwenstein Medical GmbH & Co. KG liable for any resulting property damage or personal injuries if:

- original accessories and spare parts are not used,
- the user instructions in this instructions for use are not followed,
- installation, settings, modifications, enhancements or repairs are not performed by Löwenstein Medical GmbH & Co. KG or by personnel authorized by Löwenstein Medical GmbH & Co. KG.

5 Troubleshooting

Malfunction	Display	Probable cause	Remedy
Device does not start	  do not illuminate	Battery is discharged.	Connect the power supply unit for charging or operation.
	 blinks red	Battery is discharged and impermissible input voltage.	Connect the correct power supply unit, check the power supply.
	 blinks red	Internal error	Please contact service!
Device does not charge	 blinks red and  illuminates green	- Battery temperature too high for charging - Battery is defective.	- No measure required, automatic resumption of charging when battery has cooled down. - Please contact service!
	 blinks red	Impermissible input voltage	Connect the correct power supply unit, check the power supply.
	 illuminates green/orange and/or  illuminates green	- Bacterial filter is blocked. - Leak in the suction line or in the canister - Service life of vacuum pump expired	- Replace the bacterial filter. - Check that secretion canister lid and tube connections are tight and no leakages are apparent. - Please contact service!
No flow rate	   continuous red light	- Float sphere of the overflow protection blocks the aspiration lid inlet. - Bacterial filter is blocked. - Suction tube is blocked. - Bacterial filter is not inserted.	- Check fluid level in the secretion canister, empty the canister. - Replace the bacterial filter. - Replace the suction tube. - Insert the bacterial filter.
Device switches off	 blinks red	Internal error	Please contact service!

6 Transport, storage and disposal

6.1 Decontamination prior to shipment

Before the **secretion aspirator Allegra M30** is passed on, the device must be properly processed by Löwenstein Medical GmbH & Co. KG or by specialized personnel authorized by Löwenstein Medical GmbH & Co. KG to protect subsequent users. Processing is mandatory in accordance with the German Medical Device Operator Ordinance (MPBetreibV), German Medical Devices Act (MPG) and the manufacturer's instructions.

Löwenstein Medical GmbH & Co. KG offers its partners and customers fast and professional processing as well as required testing services.

Before the **secretion aspirator Allegra M30** is shipped to Löwenstein Medical GmbH & Co. KG, it must be cleaned and disinfected. Please follow the instructions in section 4.1 resp. 4.2! Please affix the supplied "Used Medical Device" label to the shipping carton! Please notify Löwenstein Medical GmbH & Co. KG of the return in advance.

6.2 Storage

Store the **secretion aspirator Allegra M30** in accordance with the indications in the technical data (section 7)!

Charge the battery of the **secretion aspirator Allegra M30** before storing the device. This ensures that it is operational at all times.

If the **secretion aspirator Allegra M30** is not used for a longer period of time (approximately 10 months), the battery must be fully recharged again!

6.3 Disposal

- The components of the device must be disposed of properly at the end of the product's service life.
- Ensure that the disposed components are clean and carefully sorted by material.
- The housing material is labeled with a material symbol mark and is fully recyclable.
- Decontaminate the device and the accessories prior to disposal.
- According to EU Directives 2012/19/EU Waste Electrical and Electronic Equipment (WEEE) and 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS II), the device must not be disposed of in domestic waste.
- The device and the accessories may be disposed of via Löwenstein Medical GmbH & Co. KG or the service partner.
- In relation to proper disposal of potential biohazards, contact your local authority.
- Outside the EU: Follow the disposal requirements of your country!

7 Technical data

Model name	Allegra M30
Flow rate* (measuring point at suction tube nozzle)	30 l/min ± 2 l/min (high flow)
Vacuum	max. -90 kPa (~ -900 mbar; -675 mmHg) (high vacuum) Conversion factor: 1 kPa ~ 7.5 mmHg; 1 kPa ~ 10 mbar
Canister	Reusable secretion canister (1000 ml)
Suction tube	PVC, length 150 cm, vacuum stable up to -0.9 bar item number: 100713-2 inner diameter 6 mm, wall thickness 2 mm 100765 inner diameter 7 mm, wall thickness 1.5 mm
Nominal voltage of the power supply unit	In: AC 100 – 240 V ~ 50-60 Hz/1.5 A Out: DC 12 V/5.0 A
Maximum load current	5.0 A
Permissible input voltage	12 V
Power consumption at 12 V	60 W
Degree of protection pursuant to IEC 60601-1	Type BF
Risk classification pursuant to 93/42/EEC, IX	Ila
Protection class pursuant to IEC 60601-1	II
IP degree of protection	IP22

CE marking	CE0044
Sound emission	without vacuum: approximately 60 dB(A)
Ambient conditions	Transport/Storage ambient temperature: -25 °C to +70 °C relative humidity: up to 90 %, non-condensing air pressure: 700 hPa to 1060 hPa
	Operation ambient temperature: +5 °C to +40 °C relative humidity: 15 % to 90 %, non-condensing air pressure: 700 hPa to 1060 hPa
Battery	10.8 V; 5800 mAh (lithium ion battery)
Charging time if battery is empty	~ 3h
Charging time if battery is approximately 50% full	~ 2h
Power supply unit	GTM91099-6015-3.0-T2, GlobTek Inc.
Dimensions basic device (H x W x D)	213 mm x 260 mm x 124 mm
Weight (basic device)	2.8 kg
Operating time	5 minutes ON; 5 minutes OFF
Runtime in battery operation	approximately 45 minutes, depending on the runtime of the motor
Lifetime	5 years
Item number (REF)	1510215hl

**The information provided may differ depending on the altitude above sea level, the prevailing air pressure and the air temperature.*

8 EMC information

WARNING: The use of the **secretion aspirator Allegra M30** directly adjacent to or stacked with other devices should be avoided, since this could lead to impermissible operation. If the use of the **secretion aspirator Allegra M30** adjacent to or stacked with other devices is required, the **secretion aspirator Allegra M30** and the other devices should be monitored in order to verify proper operation in this arrangement!

WARNING: The use of accessories and spare parts, transformers and cables for the **secretion aspirator Allegra M30** not indicated or provided by Löwenstein Medical GmbH & Co. KG may increase the electromagnetic emissions or reduce the electromagnetic immunity of the **secretion aspirator Allegra M30**, resulting in impermissible operation. No warranty is provided for damages caused by using accessories and spare parts, transformers and cables not recommended or by improper use. Only use original accessories and spare parts!

WARNING: The use of the indicated or provided accessories and spare parts, transformers and cables for other devices than the **secretion aspirator Allegra M30** may increase the electromagnetic emissions or reduce the electromagnetic immunity. No warranty is provided for damages caused by using the indicated or provided accessories and spare parts, transformers and cables with other devices or by improper use. Only use the accessories and spare parts, transformers and cables with the **secretion aspirator Allegra M30**!

WARNING: Portable and mobile RF communication equipment (incl. peripheral devices such as antenna cables and external antennas) may influence medical electrical devices and therefore should not be used within a range of 30 cm of any part of the **secretion aspirator Allegra M30** incl. its cables. Otherwise the performance of the device may be impaired.

WARNING: The **secretion aspirator Allegra M30** may influence other devices, examinations and treatments electromagnetically. For this reason, special attention should always be paid to other devices and examinations or treatments performed in parallel so that any influence is detected as quickly as possible.

The **secretion aspirator Allegra M30** meets the requirements of IEC 60601-1-2 / EN 60601-1-2 "Electromagnetic Compatibility - Medical Electrical Equipment" without deviations or restrictions. Therefore electromagnetic interference is reduced to a minimum. Follow the indicated instructions and guidelines to sustain the basic safety and the essential features of the **secretion aspirator Allegra M30** over its entire service life.

8.1 Electromagnetic environment in which the secretion aspirator Allegra M30 may be operated

The **secretion aspirator Allegra M30** is intended for operation in the electromagnetic environment specified below, in which the RF disturbances are controlled. The customer or the user of the **secretion aspirator Allegra M30** must ensure that it is operated in such an environment.

The environments for intended operation include professional healthcare institutions and the home care environment. Special environments, such as the proximity of RF surgery or MRI resp. environments in which the intensity of the EMC disturbances is high, are excluded.

Emission limits

Conducted and radiated RF emissions	CISPR 11
Harmonic distortion	IEC 61000-3-2
Voltage fluctuations/flicker	IEC 61000-3-3

Housing

Phenomenon	Test method	Immunity test levels
Electrostatic discharge	IEC 61000-4-2	±8 kV contact ±15 kV air
Radiated RF disturbances	IEC 61000-4-3	28 V/m 80 MHz to 6 GHz 80% AM at 1 kHz
Magnetic field at rated power frequency (50 Hz)	IEC 61000-4-8	30 A/m

Power supply AC

Phenomenon	Test method	Immunity test levels
Fast transient electrical disturbances/bursts	IEC 61000-4-4	±2 kV 100 kHz repetition frequency
Impulse voltages/surges	IEC 61000-4-5	±1 kV line-to-line ±2 kV line-to-ground
Conducted RF disturbances	IEC 61000-4-6	3 V 0.15 MHz to 80 MHz 6 V in ISM bands and amateur radio bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz
Voltage dips	IEC 61000-4-11	0% U_T for ½ period At 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% U_T for 1 period 70% U_T for 25 periods one-phase at 0°
Voltage interruptions	IEC 61000-4-11	0% U_T for 250 cycles

Note: U_T is the AC mains voltage prior to the application of the test levels.

Power supply DC

Phenomenon	Test method	Immunity test levels
Fast transient electrical disturbances/bursts	IEC 61000-4-4	±2 kV 100 kHz repetition frequency
Impulse voltages/surges	IEC 61000-4-5	±1 kV line-to-line ±2 kV line-to-ground
Conducted RF disturbances	IEC 61000-4-6	3 V 0.15 MHz to 80 MHz 6 V in ISM bands and amateur radio bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz
Electrical transient conduction along supply lines	ISO 7637-2	see ISO 7637-2

Patient connection

Phenomenon	Test method	Immunity test levels
Electrostatic discharge	IEC 61000-4-2	±8 kV contact ±15 kV air
Conducted RF disturbances	IEC 61000-4-6	3 V 0.15 MHz to 80 MHz 6 V in ISM bands and amateur radio bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz

Signal input and output channels

Phenomenon	Test method	Immunity test levels
Electrostatic discharge	IEC 61000-4-2	±8 kV contact ±15 kV air
Fast transient electrical disturbances/bursts	IEC 61000-4-4	±1 kV 100 kHz repetition frequency
Impulse voltages/surges	IEC 61000-4-5	±2 kV line-to-ground
Conducted RF disturbances	IEC 61000-4-6	3 V 0.15 MHz to 80 MHz 6 V in ISM bands and amateur radio bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz

8.2 Handling of electromagnetic interactions

Although electromagnetic interference of the **secretion aspirator Allegra M30** has been reduced to a minimum, electromagnetic disturbances between the **secretion aspirator Allegra M30** and other devices cannot be excluded. Therefore you should comply with the specified requirements and instructions regarding the permissible electromagnetic environment at any case and monitor the **secretion aspirator Allegra M30** in order to ensure proper operation and to prevent adverse events for patients and users. Select another location for the **secretion aspirator Allegra M30**, if the permissible electromagnetic environment cannot be ensured or if you have observed functioning of the **secretion aspirator Allegra M30** or of another device in the vicinity that is not intended.

Since electromagnetic propagation is affected by absorption and reflection from structures, objects and people and since the field strength from fixed transmitters cannot be predicted in advance with accuracy, a site survey of the electromagnetic phenomena at the location at which the **secretion aspirator Allegra M30** shall be operated should be considered to assess the existing electromagnetic environment at this location. If unusual performance is detected, additional measures may be required, for example a modified orientation or another location of the **secretion aspirator Allegra M30**. If the essential features of the **secretion aspirator Allegra M30** are impaired by electromagnetic disturbances, a reduced or non-existent flow rate and vacuum may be expected.

8.3 List of all cables and transformers replaceable by the user

Name	Specification	Maximum length
Power supply unit incl. connecting cable	Type: GTM91099-6015-3.0-T2 Manufacturer: GlobTek Technical data: 100-240 V~, 50-60 Hz, 1.5 A (in) 12 V DC, 5 A (out)	1.20 m
Mains cable	Type: H03VVH2-F Manufacturer: HAWA Technical data: 250 V~, 2.5 A	1.80 m

9 Order information accessories

Item number	Description	PU
1510215hl	secretion aspirator Allegra M30 inclusive accessories	1
1510214hl4	Reusable secretion canister Allegra M30 (1,000 ml)	1
1510214hl5	Lid for reusable secretion canister Allegra M30	1
1510214hl6	Reusable secretion canister system Allegra M30 (complete, incl. container, lid, suction tube and fingertip)	1
1510214hl12	Sphere for overflow protection	3
1510214hl3e	Bacterial filter Allegra M30	1
1510214hl7	Universal holder Allegra M30	1
100732	Stand with 5 rollers	1
1510214hl8	Bag Allegra M30	1
1510214hl9	Rinsing bottle, reusable (250 ml)	1
1510214hl10	Tube mount Allegra M30	5
1510214hl13	Power supply unit GlobTek GTM91099-6015-3.0-T2 incl. mains cable	1
1510214hl17	Car connecting cable	1
100713-2	Suction tube with funnel and fingertip (inner diameter 6 mm, wall thickness 2 mm) (unsterile)	10
100765	Suction tube with funnel and fingertip (inner diameter 7 mm, wall thickness 1.5 mm) (unsterile)	10
214063-901	Fingertip S6	1
100733	Filing tray for roller stand	1
1510214hl11	Device holder for roller stand Allegra M30	1
1510214hl16	Hose connector S6	5

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Fabricant : Löwenstein Medical GmbH & Co. KG