



# AIRCON

## Respiratory Gas Humidifier

Homecare  
Pneumology  
Neonatology  
Anaesthesia  
INTENSIVE CARE VENTILATION  
Sleep Diagnostics  
Service  
Patient Support

# Active humidification redefined

## aircon

With the Aircon respiratory gas humidifier, a new overall concept for all patient groups and areas of use is available, which is optimally adapted for the demands of the clinical sector. The humidifier Aircon fulfills medical requirements as well as economic market expectations.

The Aircon respiratory gas humidifier combines modern technology and innovative design in a high performance unit. Thus, it fulfils the medical requirements and the economic expectations of the market.

The result is an optimal physiologically conditioned respiratory gas, which protects the mucous membranes of the ventilated patient from drying out and prevents interference with the mucocilliary clearance.

### Scope of Delivery

Humidifier Aircon, 230V
Temperature probe
Heating wire adapter (inspiratory+expiratory)
Country-specific power cord

Our breathing tube systems (single and double limb) are applicable to neonates, children and adults. We offer configuration for clinical and extra-clinical use for all common ventilators.

# Everything from one source

Aircon provides a complete respiratory humidification system in combination with its specially designed humidifier chambers and breathing tube systems. It is suitable for mechanically ventilated patients requiring individual therapy with high performance and safety.

## USER-FRIENDLY

- 3.5" TFT colour display with automatic dimmer
- Logical menu navigation with symbols and pictograms
- Treatment pause function

## SAFE

- Elaborated alarm management
- Automatic water level monitoring
- Incident and alarm protocol (data exportable to PC)

## HIGH QUALITY

- Protection Class II for clinical and extra-clinical use
- Less maintenance, no hidden costs
- Manufacturing "Made in Germany"

## SMART PERFORMANCE


- 3 function modes (IV, NIV, FREE)
- Expiratory tube: adjustable heating performance
- Individual adjustment of humidification performance
- Economical energy consumption

## EFFICIENT

- Suitable for all common ventilators
- Ready for immediate use
- Extended accessory range



Heinen + Löwenstein  
**WORLDWIDE**



p-aircon-e/1 Stand: 10/2014 Subject to technical changes without notice

