



Analogtabelle Beatmungsmodi

HAMILTON-T1	Dräger Oxylog 2000 plus / 3000 plus	Weinmann Emergency Medumat Transport
DuoPAP	PC-BIPAP	BiLevel
PCV	---	PCV
APRV	---	---
(S) CMV+	VC-CMV Autoflow	---
SIMV+	VC-SIMV Autoflow	PRVC + ASB
Spont	SpnCPAP / PS	CPAP + ASB
ASV	---	---
NIV ST	PC-BIPAP (NIV)	BiLevel (NIV)
NIV	SpnCPAP / PS (NIV)	CPAP + ASB (NIV)



Heinen + Löwenstein GmbH & Co. KG
Arzbacher Straße 80
D - 56130 Bad Ems
Fon: +49 26 03/96 00 0
Fax: +49 26 03/96 00 50
info@hul.de
www.hul.de