


Umsetzung der S3-Leitlinie NIV bei ARI



Analogtabelle Beatmungsmodi

HAMILTON C1-C3, T1, G5, S1, MR1	Dräger Evita 4 / XL	Dräger Evita V500 / V300	GE Carescape R860
P-CMV	---	PC-CMV	A/C PC
P-SIMV	---	PC-SIMV	SIMV-PC
DuoPAP	BIPAP	PC-BIPAP	BiLevel
APRV	APRV	PC-APRV	APRV
Spont	ASB	SPN-CPAP/ASB	CPAP / PSV
NIV	ASB (NIV-Menü)	SPN-CPAP/ASB (NIV-Menü)	CPAP / PSV (NIV-Menü)
NIV-ST	---	---	---
(S) CMV	IPPV, IPPVass	VC-CMV	A/C VC
SIMV	SIMV	VC-SIMV	SIMV-VC
APVcmv	IPPV-Autoflow	VC-CMV-Autoflow	A/C PRVC
APVsimv	SIMV-Autoflow	VC-SIMV-Autoflow	BiLevel VG
PC-SIMV VG	---	---	SIMV-PRVC
VS	---	SPN-CPAP/VS	VS
ASV	---	---	---



Löwenstein Medical GmbH & Co. KG
 Arzbacher Straße 80
 D - 56130 Bad Ems
 Fon: +49 26 03/96 00 0
 Fax: +49 26 03/96 00 50
 info@hul.de
 www.hul.de