



ALARMS:

- Temperature
- Pressure
- Power

INDICATIONS:

- Outlet pressure
- Working hours
- Drying



TECHNICAL SPECIFICATION

Compressor efficiency	40 l.min ⁻¹ / 3 bar
Air tank capacity	2 lit
Peak Flow at pressure 3 bar	200 Lit.min ⁻¹
Filtration of air to	5 µm
Dew point depression at 40 Lit.min ⁻¹ , 20°C	5°C under the ambient temperature
Outlet connection	DISS I 160-A (3/4"-16 UNF)
Mode of the operation	Continuous - SI
Noise level	49 dB [A]
Dimensions (w x d x h)	440 x 345 x 440
Weight	30 kg

SMART 50

Small and Compact medical compressor

Homecare
Pneumology
Neonatology
Anaesthesia
INTENSIVE CARE VENTILATION
Sleep Diagnostics
Service
Patienten Support

Löwenstein Medical
Arzbacher Straße 80
D-56130 Bad Ems
Germany

fon: + 49 2603 9600-0
fax: + 49 2603 9600-50
mail: export@hul.de
Internet: www.hul.de

