



... ingenious simplicity. Nasal or Full Face. Simply easy. It could not be any simpler. Everything about our ventilation mask JOYCEeasy is designed for easy use – from simple mask adjustments directly on the patient's face to intuitive handling. Because the mask fits so well, JOYCEeasy also encourages better therapy compliance by minimizing leaks and pressure points. The quiet, diffuse exhalation system ensures a peaceful bedroom atmosphere.

Ingenious simplicity everywhere:

- Unlimited freedom of movement provided by the ball-and-socket joint
- Quiet and diffuse exhalation system 2
- Mask can be put on and taken off with just one hand
- Color-coded hook-and-loop tabs 3
- Extremely robust plastic
- Quick-release cord (only Full Face)

JOYCEeasy next FF JOYCEeasy Full Face

... the perfect duo. The right fit for everyone.

We heard you.

And we reworked JOYCEeasy Full Face to satisfy your practical needs.

JOYCEeasy next FF now features a completely redesigned large mask cushion that has a comfortable and firm fit even at high pressures.

The headgear too remains securely in place at high pressures. Thanks to the right material mix and ergonomic design, the headgear clips can be easily loosened from the mask.

Your benefit – The right mask for everyone:

JOYCEeasy next FF = normal to obese JOYCEeasy Full Face = slender to gaunt













anti-asplyxia valve, rotating sleeve WM 15897



Endoscope adapter set JOYCEeasy next FF NV/

JOYCEeasy Full Face NV

retaining ring NV,

endoscope elbow NV WM 15864

Forehead cushion

(Nasal and Full Face)

Size S WM 26803 (Nasal) WM 26853 (Full Face) WM 25104 (next FF) Size M WM 26813 (Nasal)

WM 26863 (Full Face)

WM 25114 (next FF) Size L WM 26823 (Nasal)

WM 26873 (Full Face)

WM 25124 (next FF)

WM 26801

Mask cushion

Quick-release cord

WM 26890 (Full Face)

Ball-and-socket set, vented Easy next FF retaining ring, elbow, anti-asphyxia valve, rotating sleeve WM 15944





Product name	Size S	Size M		Size L	Ī
i roduce name	5126 5	512011		JIZE L	
JOYCEeasy Nasal	WM 25510	WM 25520		WM 25530	
JOYCEeasy Full Face	WM 25910	WM 25920		WM 25930	
JOYCEeasy Full Face NV	WM 26810	WM 26820		WM 26830	
JOYCEeasy next FF	WM 25109	WM 25119		WM 25129	
JOYCEeasy next FF NV	WM 25159	WM 25169		WM 25179	

For more information on our therapy solutions, accessories and mask systems please visit weinmann-medical.com $% \mathcal{A}_{\mathrm{res}}^{\mathrm{res}}$

Weinmann develops, produces and markets diagnostic and therapeutic solutions for sleep medicine and ventilation.

Products are tailored exactly to the needs of patients and experts. Before market launch, every product goes through a certified validation process in research and development, procurement, production and quality management.

Weinmann, a Löwenstein Group company since 2013, is among the Top 3 in many of its demanding international markets, thanks to ideas, experience and innovation. Thanks too to the decades of trust that Weinmann's partners around the world have placed in our proven products Made in Germany.

CE 0197





You Tube

Heinen + Löwenstein GmbH & Co. KG Arzbacher Straße 80 56130 Bad Ems, Germany www.hul.de Weinmann Geräte für Medizin GmbH + Co. KG Kronsaalsweg 40 22525 Hamburg, Germany T +49 40 5 47 02 - 100 F +49 40 5 47 02 - 476 E info@weinmann-medical.com www.weinmann-medical.com









JOYCEeasy

Easy to Fit. Easy to Use.

A Company in the Löwenstein Group WEINMANN medical technology