



## EasyLung™

The lung simple / Neonatal Testlunge

- Technical Specification / Technische Daten

### Technical Specification

Volume	20 mL @ PEEP = 0 mbar; $P_{insp} = 30$ mbar
Compliance	0.7 mL/mbar @ $V_t$ 20 mL
Dimensions	45 mm × 150 mm
Weight	25 g
Material	Silicone, sterilizable
Order number	0217043bg

EasyLung™ Neonatal Test lung simple

The EasyLung Neonatal is the ideal tool to test neonatal and infant ventilators.



EasyLung™ Neonatal Testlunge

Die EasyLung Neonatal ist die ideale Testlunge zur Prüfung von Neonatal-Beatmungsgeräten.

### Technische Daten

Volumen	20 mL @ PEEP = 0 mbar; $P_{insp} = 30$ mbar
Compliance	0.7 mL/mbar @ $V_t$ 20 mL
Größenangaben	45 mm × 150 mm
Gewicht	25 g
Material	Silikon, sterilisierbar
Artikelnummer	0217043bg

to the product page:

The images contain optional accessories which can be ordered in addition. Please inquire further information.

